

4

3

2

1

D

D

MINIMUM OPENING DIAMETER

 $\phi 1.23$

BODY LENGTH

2.72

CLOSED OVERALL LENGTH
2.76 $\phi 1.28$

C

C

B

B

NOTES:

1. REV1.1
2. MIN. HINGE LIFE: 25 CYCLES
3. OPENING FORCE: 18-35 LBS
4. COLOR MATCHING BATCH TO BATCH USING HUNTER L+/-1.5, A+/-1.0, B+/-1.0
5. MAX. FLASH ALLOWED: 0.050IN
6. COSMETIC EXEMPTIONS: MINIMAL SCRATCHING AND SCUFFING ALLOWED BUT MAY OCCUR FROM VARYING SHIPPING AND FRIEGHT CONDITIONS

A

A



INTERTECH

PART NAME

13 DRAM

SIZE

ASTANDARD MATERIAL
POLYPROPYLENE (PP)STANDARD WEIGHT
9.71 GMS

DIMENSIONS ARE IN INCHES

DO NOT SCALE

SCALE: 1:1

SHEET 1 OF 1

4

3

2

1