






NOTES:

- 0.020" max flash allowed
- Hinge Life cycle = 25 min.
- Opening Force = 14-24 lbs.



Part	D	H
13 Dram	1.259	2.718
19 Dram	1.482	2.718
30 Dram	1.850	2.718
60 Dram	1.850	5.325



 INTERTECH Dimensions are in Inches Dimensions are REFERENCE ONLY	Drawn By	Date 1/23/2020	Approved
	Title Dram Bottles		
	Part/Asset Number 13, 19, 30, 60	Weight	
	Material PP	REV	
Scale 1:2	Page 1/1	1	